

Labycare-OT Light LED 6 Octa Reflector Ceiling 126 LED (SKSS-OL019)



OT Light LED 6 Octa Reflector Ceiling 126 LED Manufacturer by LABYCARE

Description :

An LED 6 Octa Reflector Ceiling-Mounted OT Light with 126 LEDs is engineered for major and super-specialty surgeries, delivering ultra-high intensity, uniform, and shadow-free illumination with excellent depth perception and true color accuracy.

Specifications :

Type: Ceiling-mounted, LED 6 Octa Reflector

Light Source: 126 high-intensity LEDs (approx. 21 LEDs per reflector)

Illumination: Up to 2,20,000 Lux

Color Temperature: 4,000–4,500 K (natural white)

Light Field Diameter: 180–250 mm

Color Rendering Index (CRI): ?95

Shadow Reduction: Advanced multi-octa reflector optical system

Control: Independent intensity control via sterile handle / wall panel

Heat Emission: Cool light with minimal IR & UV radiation

Power Supply: 220–240 V AC, 50/60 Hz

Applications: Major OT, cardiac, neuro, orthopedic, trauma & super-specialty surgeries

labycare

www.labycare.com | Email: labycare@gmail.com