

Labycare-OT Light LED 5 Octa Reflector Mobile 105 LED (SKSS-OL023)



OT Light LED 5 Octa Reflector Mobile 105 LED Manufacturer by LABYCARE

Description :

The OT Light LED 5 Octa Reflector Mobile (105 LED) is a high-performance floor-mounted surgical light delivering ultra-bright, uniform, and shadow-reduced illumination. Its multi-reflector LED design and flexible mobile arm ensure precise positioning and true color rendering, ideal for minor, general, and specialty surgical procedures in OTs, emergency rooms, and clinics.

Specifications :

Specification	Details
Type	Mobile / Floor-mounted LED OT Light
LED Count	63 high-intensity LEDs
Optical Design	3 Octa reflector LED configuration
Illumination (Lux)	Approx. 1,00,000–1,60,000 Lux
Color Temperature	4,000–4,500 K (natural white)
Color Rendering Index (CRI)	≥95
Light Field Diameter	170–230 mm
Shadow Reduction	Multi-reflector octa design
Control	Intensity control on light head / panel
Arm Movement	Spring-balanced flexible arm
Mobility	Heavy-duty base with lockable castors
Heat Emission	Cool LED light, minimal IR/UV
Power Supply	220–240 V AC, 50/60 Hz
Applications	Minor OT, emergency, procedure rooms

labycare

www.labycare.com | Email: labycare@gmail.com