

## Labycare-Angle Poise Examination Light (SKSS-OL038)



Angle Poise Examination Light Manufacturer by **LABYCARE** in India.

### **Description :**

Labycare Examination Mobile Lamp Angle Poise

Examination Lamp Mobile 1 Dome Aluminium is also Known as Angle Poise Lamp

These Lamps are used in Clinics for general examination of patients.

Examination Lamp is fitted with PU Base having 5 wheels

The length of the Arm is 75 cm

The dome is made of Aluminium and looks like Shade type

Height Adjustment can be done

75 - 100 W Bulb can be used in the Examination Lamp

The Voltage of this Examination Lamp is 220 V

Knock Down Construction to Save Freight

We are Manufacturer, Exporters and Suppliers of Angle Poise Examination Light at wholesale competitive prices from India

### **Specifications :**

Labycare Examination Mobile Lamp Angle Poise

Examination Lamp Mobile 1 Dome Aluminium is also Known as Angle Poise Lamp

These Lamps are used in Clinics for general examination of patients.

Examination Lamp is fitted with PU Base having 5 wheels

The length of the Arm is 75 cm

The dome is made of Aluminium and looks like Shade type

Height Adjustment can be done

75 - 100 W Bulb can be used in the Examination Lamp

The Voltage of this Examination Lamp is 220 V

Knock Down Construction to Save Freight

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)