

Labycare-Halogen Examination Light Chrome Base (SKSS-OL041)



Halogen Examination Light Chrome Base Manufacturer by **LABYCARE** in India.

Description :

Labycare Examination Mobile Lamp Halogen

Examination Lamp Mobile 1 Dome also Known as Halogen Examination Lamp

These Lamps are used in Clinics for general examination of patients.

Examination Lamp is fitted with Sheet Metal Chrome Plated Base having 5 wheels

The length of the Arm is 71 cm

Dome is equipped with handle for easy movement.

Halogen Bulb is fitted in Halogen Examination Light

Height Adjustment can be done

The Voltage of this Examination Lamp is 220 V

Knock Down Construction to Save Freight

We are Manufacturer, Exporters and Suppliers of Halogen Examination Light Chrome Base at wholesale competitive prices from India

Specifications :

Labycare Examination Mobile Lamp Halogen

Examination Lamp Mobile 1 Dome also Known as Halogen Examination Lamp

These Lamps are used in Clinics for general examination of patients.

Examination Lamp is fitted with Sheet Metal Chrome Plated Base having 5 wheels

The length of the Arm is 71 cm

Dome is equipped with handle for easy movement.

Halogen Bulb is fitted in Halogen Examination Light

Height Adjustment can be done

The Voltage of this Examination Lamp is 220 V

Knock Down Construction to Save Freight

labycare

www.labycare.com | Email: labycare@gmail.com