

Labycare-LED 21 Examination Light (SKSS-OL047)



LED 21 Examination Light Manufacturer by **LABYCARE** in India.

Description :

This Labycare 21 led examination spot OT light specifically illuminates specific areas of a patient's body during diagnostic procedures and general check-ups. It is intended to produce its useful energy mostly in the visible range. In addition to this, our offering is necessary while doing minor surgical procedures like cosmetic surgery. This product has a hefty base and a low center of gravity fitted into it. 21 led examination spot OT light is equipped with five anti-static swivel castors mounted on a base in the form of a star. It is prepared using polyvinyl chloride material and modern technology.

We are Manufacturer, Exporters and Suppliers of LED 21 Examination Light at wholesale competitive prices from India.

Specifications :

Suitable for Minor Surgeries, Examination, Stitches etc, Ease of Use.

* Available in Wall Mounted and Mobile Stand.

Max Light Intensity (@1M) : 40000 lux

Dimming Range : Fixed

Field Diameter (Spot) : 150 mm

Diameter of Dome : 16 mm

No. of LEDs : 21

Focusing : Fixed

LED Life Time : 50,000 HRS

CRI : >92

Colour Temperature : 4,300K

Depth of Illumination : 80 cm

Power Consumption : 12 Watt

Power Supply : 220V/50 Hz AC

Temp. Rise At Surgeon's Head : < 1°C.

labycare

www.labycare.com | Email: labycare@gmail.com