

## Labycare-Halogen Ceiling OT Light Single Dome (SKSS-OL049)



Halogen Ceiling OT Light Single Dome Manufacturer by **LABYCARE** in India.

### **Description :**

Labycare Halogen Ceiling Surgical Lamp Single Dome

Dome Dia : 500 mm

Handle is Sterilizable

Lux intensity : 50000  $\pm$  10%

Halogen Bulb of 24V/150W

Heat Absorbing Light

All fittings are made of SS.

Colour correcting filter glass

In case of Emergency, Reserve bulb is there

Light adjustment can be increase or decrease through dimmer.

Knock Down Construction to Save Freight

We are Manufacturer, Exporters and Suppliers of Ceiling OT Light Single Dome at wholesale competitive prices from India.

### **Specifications :**

Labycare Halogen Ceiling Surgical Lamp Single Dome

Dome Dia : 500 mm

Handle is Sterilizable

Lux intensity :  $50000 \pm 10\%$

Halogen Bulb of 24V/150W

Heat Absorbing Light

All fittings are made of SS.

Colour correcting filter glass

In case of Emergency, Reserve bulb is there

Light adjustment can be increase or decrease through dimmer.

Knock Down Construction to Save Freight

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)