

# Labycare-Hand Set Heat Sealer (SKSSHSHS947)



Hand Set Heat Sealer

## **Description :**

Hand Set Heat Sealer

<b><u>Specifications :</u></b>	Hand Set Heat Sealer	Hand Set Heat Sealer
<b><u>TECHNICAL</u></b>	Material	Plastic
<b><u>SPECIFICATIONS:</u></b>	Plug Type	US Stand
	Suitable For	PE/PP E
	Power	300 W
	Size	8 Inch
	Sealing Length	20 cm
	Weight	1300 g
	Voltage	110 V,50
	Color	Blue
	Heat Time	0.2 -1.3s
	Dimensions	L12.99x
	Max Sealing Thickness	0.2 mm

labycare

www.labycare.com | Email: labycare@gmail.com