

Labycare-Infant Feeding Tube with Graduation (SKSS-MD030)



Infant Feeding Tube with Graduation Manufacturer by **LABYCARE** in India.

Description :

Labycare Infant Feeding Tube is made from non-toxic, medical grade PVC.

Feeding Tube is used for neonates and paediatric feeding.

Feeding Tube has smooth coned distal end for atraumatic intubation, with two lateral eyes for efficient drainage.

Feeding Tube has markings at 20 cm and radio opaque line for accurate placement and X-ray visualization.

Feeding Tube has proximal end with luer mount for connection with feeding device.

Marking from the tip 5 cm to 40 cm.

We are Manufacturer, Exporters and Suppliers of Infant Feeding Tube with Graduation at wholesale competitive prices from India.

Specifications :

Labycare Infant Feeding Tube is made from non-toxic, medical grade PVC.

Feeding Tube is used for neonates and paediatric feeding.

Feeding Tube has smooth coned distal end for atraumatic intubation, with two lateral eyes for efficient drainage.

Feeding Tube has markings at 20 cm and radio opaque line for accurate placement and X-ray visualization.

Feeding Tube has proximal end with luer mount for connection with feeding device.

Marking from the tip 5 cm to 40 cm.

labycare

www.labycare.com | Email: labycare@gmail.com