

Labycare-Ryle's Tube Elastomer Coated (SKSS-MD033)



Ryle's Tube Elastomer Coated Manufacturer by **LABYCARE** in India.

Description :

Labycare Ryle's Tube is made from non-toxic, medical grade PVC and is Elastomer Coated.

Ryle's Tube is used for nutritional feeding purposes or aspiration of intestinal secretions for enhanced patient safety.

Ryle's Tube has coned distal end with steel balls sealed into the tube to facilitate easy insertion.

Ryle's Tube has four lateral eyes for efficient drainage.

Ryle's Tube has proximal end universal funnel connector.

Ryle's Tube has markings at 50, 60 and 70 cm and radio opaque line for accurate placement.

Effective Length is 119 cm

Total Length is 122 cm

We are Manufacturer, Exporters and Suppliers of Ryle's Tube Elastomer Coated at wholesale competitive prices from India.

Specifications :

Labycare Ryle's Tube is made from non-toxic, medical grade PVC and is Elastomer Coated.

Ryle's Tube is used for nutritional feeding purposes or aspiration of intestinal secretions for enhanced patient safety.

Ryle's Tube has coned distal end with steel balls sealed into the tube to facilitate easy insertion.

Ryle's Tube has four lateral eyes for efficient drainage.

Ryle's Tube has proximal end universal funnel connector.

Ryle's Tube has markings at 50, 60 and 70 cm and radio opaque line for accurate placement.

Effective Length is 119 cm

Total Length is 122 cm

labycare

www.labycare.com | Email: labycare@gmail.com