

Labycare-Salem Sump Tube (SKSS-MD034)



Salem Sump Tube Manufacturer by **LABYCARE** in India.

Description :

Labycare Salem Sump Tube is made from non-toxic medical grade PVC

Salem Sump Tube is used to prevent the accumulation of gastric secretions and air

Double lumen tube throughout the length provides airvent when the stomach has been evacuated of accumulated liquid, air or small particulate matter

Salem Sump Tube has coned distal end with ten lateral eyes for faster drainage

Salem Sump Tube is fitted with universal funnel connector along with double taper connector at proximal end for easy extension

Salem Sump Tube has markings at 45, 55, 65 and 75 cm and radio opaque line for accurate placement

Effective Length is 119 cm

Total Length is 122 cm

We are Manufacturer, Exporters and Suppliers of Salem Sump Tube at wholesale competitive prices from India.

Specifications :

Labycare Salem Sump Tube is made from non-toxic medical grade PVC

Salem Sump Tube is used to prevent the accumulation of gastric secretions and air

Double lumen tube throughout the length provides airvent when the stomach has been evacuated of

accumulated liquid, air or small particulate matter

Salem Sump Tube has coned distal end with ten lateral eyes for faster drainage

Salem Sump Tube is fitted with universal funnel connector along with double taper connector at proximal end for easy extension

Salem Sump Tube has markings at 45, 55, 65 and 75 cm and radio opaque line for accurate placement

Effective Length is 119 cm

Total Length is 122 cm

labycare

www.labycare.com | Email: labycare@gmail.com