

Labycare-Suction Catheter Thumb Control (SKSS-MD039)



Suction Catheter Thumb Control Manufacturer by **LABYCARE** in India.

Description :

Labycare Suction Catheter is made from non-toxic, non-irritant medical grade PVC.

Suction Catheter is used for removal of secretions from the patient's mouth, trachea and bronchial tube.

Lateral eye at distal end for effective suctioning.

Proximal end of tube provided with universal connector for easy connection with suction source.

Available with graduation.

Effective Length is 480mm

Total Length is 530mm

We are Manufacturer, Exporters and Suppliers of Suction Catheter Thumb Control at wholesale competitive prices from India.

Specifications :

Labycare Suction Catheter is made from non-toxic, non-irritant medical grade PVC.

Suction Catheter is used for removal of secretions from the patient's mouth, trachea and bronchial tube.

Lateral eye at distal end for effective suctioning.

Proximal end of tube provided with universal connector for easy connection with suction source.

Available with graduation.

Effective Length is 480mm

Total Length is 530mm

labycare

www.labycare.com | Email: labycare@gmail.com