

Labycare-Urethral Catheter (SKSS-MD043)



Urethral Catheter Manufacturer by **LABYCARE** in India.

Description :

Labycare Urethral Catheter is made from non-toxic, non-irritant, medical grade PVC.

Urethral Catheter is used for one time urine drainage for patients.

Radio opaque line provided through out the catheter for X-ray visualization.

Distal end with lateral eyes and proximal end funnel connector, for maximum drainage and easy connection to the urine bag.

We are Manufacturer, Exporters and Suppliers of Urethral Catheter at wholesale competitive prices from India.

Specifications :

Effective Length (Male) is 370mm

Effective Length (Female) is 170mm

Total Length (Male) is 400mm

Total Length Female is 200mm

labycare

www.labycare.com | Email: labycare@gmail.com