

Labycare-Colony Counter (SKSS870)



Our digital colony counter is a feature packed ft easy to use, tins is an indispensable bench top tool for the busy microbiologist.It is designed for rapid and accurate counting of bacterial and mould colonies.

Description :

Digital colony counter is designed for quick and accurate counting of bacterial and mould colonies in petri dishes. Our digital colony counter is a feature packed ft easy to use, tins is an indispensable bench top tool for the busy microbiologist. It is designed for rapid and accurate counting of bacterial and mould colonies.

Specifications :

Technical Specifications

Display	:	4 Digit, 9999 Max. count
Dish Size	:	110 mm
Magnification	:	1.7
Dimension	:	L 274 X B 320 X H 167 mm (Approx)
Weight :		3 Kgs. (Approx)

Power : 230V – 10% AC, 50 Hz, 40W

Accessories : Marking Pen – 1 No.

Magnifier Lens – 1 No.

labycare

www.labycare.com | Email: labycare@gmail.com