

Labycare-Thoracic Trocar Catheter (SKSS-MD044)



Thoracic Trocar Catheter Manufacturer by **LABYCARE** in India.

Description :

Labycare Thoracic Trocar Catheter is made from non-toxic, medical grade PVC.

Thoracic Trocar Catheter is used for rapid non-operative plueral and chest drainage.

Catheter is snuggly fixed over the trocar, and provided with large smooth eyes for efficient drainage.

Catheter has smooth taper shape tip.

Catheter has Markings at 5, 10, 15 and 20cm from the distal end to measure depth of placement.

Proximal end fitted with taper shape connecter for easy connection to drainage bottle/bag.

We are Manufacturer, Exporters and Suppliers of Thoracic Trocar Catheter at wholesale competitive prices from India.

Specifications :

Labycare Thoracic Trocar Catheter is made from non-toxic, medical grade PVC.

Thoracic Trocar Catheter is used for rapid non-operative plueral and chest drainage.

Catheter is snuggly fixed over the trocar, and provided with large smooth eyes for efficient drainage.

Catheter has smooth taper shape tip.

Catheter has Markings at 5, 10, 15 and 20cm from the distal end to measure depth of placement.

Proximal end fitted with taper shape connecter for easy connection to drainage bottle/bag.

labycare

www.labycare.com | Email: labycare@gmail.com