

Labycare-Thoracic / Chest Drainage Catheter Straight (SKSS-MD045)



Thoracic / Chest Drainage Catheter Straight Manufacturer by **LABYCARE** in India.

Description :

Labycare Chest Drainage Catheter is used for post-operative drainage after cardio thoracic surgery.

Thoracic Drainage Catheter is made from non-toxic, medical grade PVC.

Chest Drainage Catheter is Smooth and round open distal end with large eyes for maximum drainage.

Thoracic Drainage Catheter is marked at every 2 cm from the last eye to ascertain the depth of placement.

Catheter is Proximal end fitted with taper tongue connector for better forceps grip and smooth penetration.

We are Manufacturer, Exporters and Suppliers of Thoracic / Chest Drainage Catheter Straight at wholesale competitive prices from India.

Specifications :

Labycare Chest Drainage Catheter is used for post-operative drainage after cardio thoracic surgery.

Thoracic Drainage Catheter is made from non-toxic, medical grade PVC.

Chest Drainage Catheter is Smooth and round open distal end with large eyes for maximum drainage.

Thoracic Drainage Catheter is marked at every 2 cm from the last eye to ascertain the depth of placement.

Catheter is Proximal end fitted with taper tongue connector for better forceps grip and smooth penetration.

labycare

www.labycare.com | Email: labycare@gmail.com