

Labycare-Manual Vacuum Aspiration Kit (SKSS-MD050)



Manual Vacuum Aspiration Kit Manufacturer by **LABYCARE** in India.

Description :

Labycare Manual Vacuum Aspiration Kit is used for Manual Aspiration procedure (Medical termination of pregnancy).

It is made from medical grade polypropylene.

The Kit consists of MTP Syringe, Karman Cannula, lubricating oil.

The MTP Syringe provides high vacuum suction during abortion procedure.

Smooth cannula with rounded tip and smooth eyes for atraumatic cannulation.

Double valve syringe.

Karman Cannula size 4mm – 12mm

We are Manufacturer, Exporters and Suppliers of Manual Vacuum Aspiration Kit at wholesale competitive prices from India.

Specifications :

Labycare Manual Vacuum Aspiration Kit is used for Manual Aspiration procedure (Medical termination of pregnancy).

It is made from medical grade polypropylene.

The Kit consists of MTP Syringe, Karman Cannula, lubricating oil.

The MTP Syringe provides high vacuum suction during abortion procedure.

Smooth cannula with rounded tip and smooth eyes for atraumatic cannulation.

Double valve syringe.

Karman Cannula size 4mm – 12mm

labycare

www.labycare.com | Email: labycare@gmail.com