

Labycare-Abdominal Drainage Kit (SKSS-MD052)



Abdominal Drainage Kit Manufacturer by **LABYCARE** in India.

Description :

Labycare Abdominal Drainage Kit is used for abdominal drainage.

It is made from non-toxic, medical grade PVC.

Drainage Kit consists of smooth and round distal end with large atraumatic eyes for efficient drainage.

Radio-opaque line provided on the catheter for x-ray visualization.

Drainage Kit is marked at every 2 cm from the cast eye to ascertain the depth of placement.

Proximal end is provided with universal taper connector for easy connection to drainage bag.

It has 2000 ml capacity bag with graduated scale.

We are Manufacturer, Exporters and Suppliers of Abdominal Drainage Kit at wholesale competitive prices from India.

Specifications :

- Used for abdominal drainage
- Made from non-toxic, medical grade PVC
- Smooth and round distal end with large atraumatic eyes for efficient drainage
- Radio-opaque line provided on the catheter for x-ray visualization
- Marked at every 2cm from the cast eye to ascertain the depth of placement
- Proximal end provided with universal taper connector for easy connection to drainage bag
- 2000 ml capacity bag with graduated scale

labycare

www.labycare.com | Email: labycare@gmail.com