

Labycare-Infusion Set Standard (SKSS-MD055)



Infusion Set Standard Manufacturer by **LABYCARE** in India.

Description :

Labycare Infusion Set is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Infusion Set has precision moulded drip chamber with sharp spike for smooth penetration.

Infusion Set has clear transparent 14mm chamber to see fluid flow.

It is provided with roller clamp for controlled flow of fluid.

Infusion Set has self sealing Latex tube / Moulded tube for additional medication.

Infusion Set is sharp, smooth vein needle (21 G) for non traumatic cannulation

Gauge Colour is Black

Infusion Set is sterile by ETO gas.

The set is for single use, non-toxic, Non Pyrogenic and sterile.

Infusion Set is individually packed in plastic / paper pouch.

We are Manufacturer, Exporters and Suppliers of Infusion Set Standard at wholesale competitive prices from India.

Specifications :

Labycare Infusion Set is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Infusion Set has precision moulded drip chamber with sharp spike for smooth penetration.

Infusion Set has clear transparent 14mm chamber to see fluid flow.

It is provided with roller clamp for controlled flow of fluid.

Infusion Set has self sealing Latex tube / Moulded tube for additional medication.

Infusion Set is sharp, smooth vein needle (21 G) for non traumatic cannulation

Gauge Colour is Black

Infusion Set is sterile by ETO gas.

The set is for single use, non-toxic, Non Pyrogenic and sterile.

Infusion Set is individually packed in plastic / paper pouch.

labycare

www.labycare.com | Email: labycare@gmail.com