

Labycare-Infusion Set Deluxe (SKSS-MD056)



Infusion Set Deluxe Manufacturer by **LABYCARE** in India.

Description :

Labycare Infusion Set Deluxe is Smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Infusion Set is fitted with precision moulded drip chamber with sharp spike for smooth penetration.

15mm Clear transparent chamber allows to see fluid flow.

It is provided with roller clamp for controlled flow of fluid.

It has self sealing Latex tube / Moulded tube for additional medication.

Vein Needle (21G) is sharp and smooth for non traumatic cannulation.

Gauge Colour is Green

Infusion Set is sterile by ETO gas.

Infusion Set is for single use, non-toxic, Non Pyrogenic and sterile.

Infusion Set is individually packed in plastic / paper pouch.

We are Manufacturer, Exporters and Suppliers of Infusion Set Deluxe at wholesale competitive prices from India.

Specifications :

Labycare Infusion Set Deluxe is Smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Infusion Set is fitted with precision moulded drip chamber with sharp spike for smooth penetration.

15mm Clear transparent chamber allows to see fluid flow.

It is provided with roller clamp for controlled flow of fluid.

It has self sealing Latex tube / Moulded tube for additional medication.

Vein Needle (21G0 is sharp and smooth for non traumatic cannulation.

Gauge Colour is Green

Infusion Set is sterile by ETO gas.

Infusion Set is for single use, non-toxic, Non Pyrogenic and sterile.

Infusion Set is individually packed in plastic / paper pouch.

labycare

www.labycare.com | Email: labycare@gmail.com