

Labycare-Infusion Set with Builtin Airvent (SKSS-MD057)



Infusion Set with Builtin Airvent Manufacturer by **LABYCARE** in India.

Description :

Labycare Infusion Set with Builtin Airvent is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Infusion Set has Precision moulded drip chamber with builtin Airvent.

Infusion Set is provided with roller clamp for controlled flow of fluid.

Self sealing Latex tube / Moulded tube for additional medication.

Vein Needle (21 G) is Sharp, smooth for non traumatic cannulation.

Infusion Set is sterile by ETO gas.

Infusion Set is for single use, non-toxic, Non Pyrogenic and sterile.

Infusion Set is individually packed in plastic / paper pouch.

We are Manufacturer, Exporters and Suppliers of Infusion Set with Builtin Airvent at wholesale competitive prices from India.

Specifications :

Labycare Infusion Set with Builtin Airvent is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Infusion Set has Precision moulded drip chamber with builtin Airvent.

Infusion Set is provided with roller clamp for controlled flow of fluid.

Self sealing Latex tube / Moulded tube for additional medication.

Vein Needle (21 G) is Sharp, smooth for non traumatic cannulation.

Infusion Set is sterile by ETO gas.

Infusion Set is for single use, non-toxic, Non Pyrogenic and sterile.

Infusion Set is individually packed in plastic / paper pouch.

labycare

www.labycare.com | Email: labycare@gmail.com