

Labycare-Blood Transfusion Set (SKSS-MD058)



Blood Transfusion Set Manufacturer by **LABYCARE** in India.

Description :

Labycare Blood Transfusion Set is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Blood Transfusion Set is fitted with Precision moulded drip chamber and sharp spike for smooth penetration.

Blood Transfusion Set has clear transparent chamber allows to see fluid flow.

Blood Transfusion Set is provided with roller clamp for controlled flow of fluid.

Self sealing Latex tube / Moulded tube for additional medication.

Vein needle (18 G) is Sharp, smooth for non traumatic cannulation.

Sterile by ETO gas.

The set is for single use, non-toxic, Non Pyrogenic and sterile.

We are Manufacturer, Exporters and Suppliers of Blood Transfusion Set at wholesale competitive prices from India.

Specifications :

Labycare Blood Transfusion Set is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Blood Transfusion Set is fitted with Precision moulded drip chamber and sharp spike for smooth penetration.

Blood Transfusion Set has clear transparent chamber allows to see fluid flow.

Blood Transfusion Set is provided with roller clamp for controlled flow of fluid.

Self sealing Latex tube / Moulded tube for additional medication.

Vein needle (18 G) is Sharp, smooth for non traumatic cannulation.

Sterile by ETO gas.

The set is for single use, non-toxic, Non Pyrogenic and sterile.

We are Manufacturer, Exporters and Suppliers of Blood Transfusion Set at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com