

Labycare-Microdrip Infusion Set (SKSS-MD059)



Microdrip Infusion Set Manufacturer by **LABYCARE** in India.

Description :

Labycare Microdrip Infusion Set is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Microdrip Infusion Set is fitted with Precision moulded drip chamber with sharp spike for smooth penetration.

Stainless Steel Cannula is provided in the chamber for controlled flow of the fluid.

Clear transparent chamber allows to see fluid flow.

Provided roller clamp for controlled flow of fluid.

Self sealing Latex tube for additional medication.

Vein needle (23 G) is sharp, smooth for non traumatic cannulation.

Microdrip Infusion Set is Sterile by ETO gas.

Microdrip Infusion Set is for single use, non-toxic, Non Pyrogenic and sterile.

We are Manufacturer, Exporters and Suppliers of Microdrip Infusion Set at wholesale competitive prices from India.

Specifications :

Labycare Microdrip Infusion Set is smooth, flexible, kink resistant tubing which is made up of non toxic medical grade P.V.C. granules.

Microdrip Infusion Set is fitted with Precision moulded drip chamber with sharp spike for smooth

penetration.

Stainless Steel Cannula is provided in the chamber for controlled flow of the fluid.

Clear transparent chamber allows to see fluid flow.

Provided roller clamp for controlled flow of fluid.

Self sealing Latex tube for additional medication.

Vein needle (23 G) is sharp, smooth for non traumatic cannulation.

Microdrip Infusion Set is Sterile by ETO gas.

Microdrip Infusion Set is for single use, non-toxic, Non Pyrogenic and sterile.

We are Manufacturer, Exporters and Suppliers of Microdrip Infusion Set at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com