

## Labycare-Scalp Vein Set (SKSS-MD060)



Scalp Vein Set Manufacturer by **LABYCARE** in India.

### **Description :**

Labycare Scalp Vein Set is also known as Butterfly set used for long term infusion

Scalp vein set is designed to provide rapid venous access with greater patient comfort during infusion.

It has short beveled siliconised needle for non traumatic cannulation.

Thin wall needle provides better flow rate per gauge, as in the same cannula more circumference is available for better flow of fluid.

Scalp Vein Set is sterile by ETO gas.

Scalp Vein Set is for single use, non-toxic, Non Pyrogenic and sterile.

Scalp Vein Set is individually packed in plastic / blister pack.

We are Manufacturer, Exporters and Suppliers of Scalp Vein Set at wholesale competitive prices from India.

### **Specifications :**

Labycare Scalp Vein Set is also known as Butterfly set used for long term infusion

Scalp vein set is designed to provide rapid venous access with greater patient comfort during infusion.

It has short beveled siliconised needle for non traumatic cannulation.

Thin wall needle provides better flow rate per gauge, as in the same cannula more circumference is available for better flow of fluid.

Scalp Vein Set is sterile by ETO gas.

Scalp Vein Set is for single use, non-toxic, Non Pyrogenic and sterile.

Scalp Vein Set is individually packed in plastic / blister pack.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)