

Labycare-Burette Set (Measured Volume Paediatric Set) (SKSS-MD061)



Burette Set (Measured Volume Paediatric Set) Manufacturer by **LABYCARE** in India.

Description :

Labycare Burette Set is also known as Measured Volume Paediatric Set

It contains Cylinder type Measure Volume chambers with 1ml increments and 10ml overflow limit.

It is provided with rubber stopper to stop fluid insertion automatically when chamber is fully filled.

Roller type flow controller provides accurate flow control.

Special plug provided at the top of the cylinder to inject the extra drug.

It is provided with spare airway.

Burette Set is Sterile by ETO gas.

Burette Set is Non toxic, Non Pyrogenic.

We are Manufacturer, Exporters and Suppliers of Burette Set (Measured Volume Paediatric Set) at wholesale competitive prices from India.

Specifications :

Labycare Burette Set is also known as Measured Volume Paediatric Set

It contains Cylinder type Measure Volume chambers with 1ml increments and 10ml overflow limit.

It is provided with rubber stopper to stop fluid insertion automatically when chamber is fully filled.

Roller type flow controller provides accurate flow control.

Special plug provided at the top of the cylinder to inject the extra drug.

It is provided with spare airway.

Burette Set is Sterile by ETO gas.

Burette Set is Non toxic, Non Pyrogenic.

We are Manufacturer, Exporters and Suppliers of Burette Set (Measured Volume Paediatric Set) at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com