

Labycare-Intravenous Cannula (SKSS-MD062)



Intravenous Cannula Manufacturer by **LABYCARE** in India.

Description :

Labycare Intravenous Cannula is also known as IV Cannula or I V Cannula

Intravenous Cannula is made of Stainless Steel

IV Cannula has double back cut ground stainless steel needle for smooth and painless vein puncture.

I V Cannula is siliconised.

PTFE Catheter: Internationally accepted excellent quality Teflon (Poly Tetra Fluoro Ethylene). The tip of the catheter is tapered on automatic machines for easy insertion into the vein.

Non - Return Valve: Side injection valve for easy access of intermittent administration of fluid with non-return valve which prevents the back flow.

Peel off Paper: A medical grade grid paper is used for blister packing which provides the effect of easy & smooth peeling

Slanted & flexible wings for safe & simple fixation.

Sterile by ETO gas.

Non toxic. Non Pyrogenic.

We are Manufacturer, Exporters and Suppliers of Intravenous Cannula at wholesale competitive prices from India.

Specifications :

Labycare Intravenous Cannula is also known as IV Cannula or I V Cannula

Intravenous Cannula is made of Stainless Steel

IV Cannula has double back cut ground stainless steel needle for smooth and painless vein puncture.

I V Cannula is siliconised.

PTFE Catheter: Internationally accepted excellent quality Teflon (Poly Tetra Fluoro Ethylene). The tip of the catheter is tapered on automatic machines for easy insertion into the vein.

Non - Return Valve: Side injection valve for easy access of intermittent administration of fluid with non-return valve which prevents the back flow.

Peel off Paper: A medical grade grid paper is used for blister packing which provides the effect of easy & smooth peeling

Slanted & flexible wings for safe & simple fixation.

Sterile by ETO gas.

Non toxic. Non Pyrogenic.

We are Manufacturer, Exporters and Suppliers of Intravenous Cannula at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com