

Labycare-Dispovan Insulin Syringe (SKSS-MD068)



Dispovan Insulin Syringe Manufacturer by **LABYCARE** in India.

Description :

Labycare Barrels are made of nontoxic, clear medical-grade polypropylene compatible with any medication.

Plungers, Needle Covers and Bottom Caps (Multipack) are made of nontoxic, medical-grade polypropylene compatible with any medication.

Needle is integrated with the barrel.

Dispovan Unolok Gaskets, Luer Lock

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

Luer Lock nozzle present in unolok syringes ensures secure fitment of the needle.

Colours Available : Red & Orange.

We are Exporters and Suppliers of Dispovan Insulin Syringe at wholesale competitive prices from India

Specifications :

Labycare Barrels are made of nontoxic, clear medical-grade polypropylene compatible with any medication.

Plungers, Needle Covers and Bottom Caps (Multipack) are made of nontoxic, medical-grade polypropylene compatible with any medication.

Needle is integrated with the barrel.

Dispovan Unolok Gaskets, Luer Lock

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

Luer Lock nozzle present in unolok syringes ensures secure fitment of the needle.

Colours Available : Red & Orange.

We are Exporters and Suppliers of Dispovan Insulin Syringe at wholesale competitive prices from India

labycare

www.labycare.com | Email: labycare@gmail.com