

Labycare-Dispovan Tuberculin Syringe (SKSS-MD069)



Dispovan Tuberculin Syringe Manufacturer by **LABYCARE** in India.

Description :

Labycare 1ml Tuberculin syringes are long type syringe with 26 Gx 1/2 needle on any other needle size as per the requirement.

These syringes are used to give BCG Vaccination.

Barrels & Plungers are made of nontoxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application.

The nozzle of the Barrels is luermount type having a 6% luer taper.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

We are Exporters and Suppliers of Dispovan Tuberculin Syringe at wholesale competitive prices from India.

Specifications :

Labycare 1ml Tuberculin syringes are long type syringe with 26 Gx 1/2 needle on any other needle size as per the requirement.

These syringes are used to give BCG Vaccination.

Barrels & Plungers are made of nontoxic, medical-grade polypropylene compatible with the medication for short term contact drug delivery application.

The nozzle of the Barrels is luermount type having a 6% luer taper.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short

term contact drug delivery application.

We are Exporters and Suppliers of Dispovan Tuberculin Syringe at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com