

Labycare-Dispovan Syringe without Needle (SKSS-MD072)



Dispovan Syringe without Needle Manufacturer by **LABYCARE** in India.

Description :

Labycare Luer mount Single use syringes without Needle.

Barrels are made of nontoxic, medical-grade polypropylene compatible with any medication.

The tip of the barrel is luer mount type having a 6% luer taper.

Plungers are made of nontoxic, medical-grade polypropylene compatible with any medication.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

We are Exporters and Suppliers of Dispovan Syringe without Needle at wholesale competitive prices from India.

Specifications :

Labycare Luer mount Single use syringes without Needle.

Barrels are made of nontoxic, medical-grade polypropylene compatible with any medication.

The tip of the barrel is luer mount type having a 6% luer taper.

Plungers are made of nontoxic, medical-grade polypropylene compatible with any medication.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

We are Exporters and Suppliers of Dispovan Syringe without Needle at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com