

Labycare-Unolock Syringe without Needle (SKSS-MD073)



Unolock Syringe without Needle Manufacturer by **LABYCARE** in India.

Description :

Labycare Luer lock Single use syringes without needle.

Secure fitment of needle due to presence of luer lock.

High clarity plunger and barrel.

Better plunger retention due to locking at the base.

Barrels are made of nontoxic, medical-grade polypropylene compatible with any medication.

The tip of the barrel is luer lock having a 6% luer taper.

Plungers are made of nontoxic, medical-grade polypropylene compatible with any medication.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

Luer Lock nozzle present in UNOLOK syringes ensures secure fitment of the needle.

We are Exporters and Suppliers of Unolock Syringe without Needle at wholesale competitive prices from India.

Specifications :

Labycare Luer lock Single use syringes without needle.

Secure fitment of needle due to presence of luer lock.

High clarity plunger and barrel.

Better plunger retention due to locking at the base.

Barrels are made of nontoxic, medical-grade polypropylene compatible with any medication.

The tip of the barrel is luer lock having a 6% luer taper.

Plungers are made of nontoxic, medical-grade polypropylene compatible with any medication.

Gaskets are made of natural rubber which is chemically inert and compatible with medication for short term contact drug delivery application.

Luer Lock nozzle present in UNOLOK syringes ensures secure fitment of the needle.

We are Exporters and Suppliers of Unolock Syringe without Needle at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com