

Labycare-Closed Wound Suction Unit (SKSS-MD078)



Closed Wound Suction Unit Manufacturer by **LABYCARE** in India.

Description :

Labycare Closed wound drainage device under negative pressure.

Closed Wound Suction Unit is provided with the option to use one or two catheters post-operatively.

Connecting tube is kink resistant and has additional strength to withstand the suction.

Multiperforated catheters provided with radio opaque line and satin smooth perforation/ eyes for trauma free performance.

Available in different sizes with matching size curved needle to meet moderate to heavy drainage needs.

Graduated bellow allows user to measure the drainage volume.

Capacity : 800ml

Flexible bellow chamber is easy to be depressed by one person single handedly for activating the suction of bellow unit.

Sterile, individually packed in peelable blister pack.

We are Manufacturer, Exporters and Suppliers of Closed Wound Suction Unit at wholesale competitive prices from India.

Specifications :

Labycare Closed wound drainage device under negative pressure.

Closed Wound Suction Unit is provided with the option to use one or two catheters post-operatively.

Connecting tube is kink resistant and has additional strength to withstand the suction.

Multiperforated catheters provided with radio opaque line and satin smooth perforation/ eyes for trauma free performance.

Available in different sizes with matching size curved needle to meet moderate to heavy drainage needs.

Graduated bellow allows user to measure the drainage volume.

Capacity : 800ml

Flexible bellow chamber is easy to be depressed by one person single handedly for activating the suction of bellow unit.

Sterile, individually packed in peelable blister pack.

We are Manufacturer, Exporters and Suppliers of Closed Wound Suction Unit at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com