

Labycare-Bougies Vented (SKSS-MD093)



Bougies Vented Manufacturer by **LABYCARE** in India.

Description :

Labycare Bougies Vented provides rigidity to endotracheal tubes to facilitate endotracheal intubation

Bougies has ultra smooth surface for low friction

Bougies has atraumatic angled tip for ease of intubation which minimise risk of patient trauma

Bougies has marking along the length to facilitate identification of tube position

Bougies has vented lumen for oxygen feed.

We are Manufacturer, Exporters and Suppliers of Bougies Vented at wholesale competitive prices from India.

Specifications :

Labycare Bougies Vented provides rigidity to endotracheal tubes to facilitate endotracheal intubation

Bougies has ultra smooth surface for low friction

Bougies has atraumatic angled tip for ease of intubation which minimise risk of patient trauma

Bougies has marking along the length to facilitate identification of tube position

Bougies has vented lumen for oxygen feed.

We are Manufacturer, Exporters and Suppliers of Bougies Vented at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com