

Labycare-Digital Audiometer Nano (SKSS-DP001)



Digital Audiometer Nano Manufacturer by **LABYCARE** in India.

Description :

A Digital Audiometer Nano is a compact, portable hearing assessment device used to evaluate hearing thresholds accurately. Designed for clinics, hospitals, and screening programs, it offers reliable pure-tone audiometry with user-friendly controls, making it suitable for routine hearing tests and early detection of hearing loss.

Specifications :

Input Voltage	AC 230 Volts - 50 Hz / DC 12 V
Air Conduction	25 Hz to 8 KHz (Maximum 120
Speech	70 dB
Bone Conduction	250 Hz to 6 KHz (Maximum 70
Speech	70 dB
Noise	Wide, Narrow & Speech Band.
Attenuator (Masking)	Steps of 5 dB
Masking Range	0 db to 100 dB

labycare

www.labycare.com | Email: labycare@gmail.com