

Labycare-ENT Diagnostic Set Basic (SKSS-DP0031)



ENT Diagnostic Set Basic Manufacturer by **LABYCARE** in India.

Description :

Labycare ENT Set is known as Diagnostic Set. It is used for physical examination of Ear, Nose and Throat.

We are Manufacturer, Exporters and Suppliers of ENT Diagnostic Set Basic at wholesale competitive prices from India.

Specifications :

Labycare ENT Set Basic includes :

- 1 double lens otoscope head with 3 reusable speculas
- 1 Ophthalmoscope head with lens selection wheel (+0 to +20 and -0 to -30) diopters single large circle
- 1 Bent Arm illuminator & 1 Nasal Speculum
- 1 Battery handle with rheostat button
- 2 laryngeal mirrors

labycare

www.labycare.com | Email: labycare@gmail.com