

Labycare-Glucometer Accu chek Active (SKSS-DP0039)



Glucometer Accu chek Active Manufacturer by **LABYCARE** in India.

Description :

Labycare Accu-Chek Active Blood Glucose Meter packs has many powerful features.

No coding reduces handling steps, you can see test result averages up to 90 days, get alerts if your test strips are expired or if you've under dosed, re-dose within 10 seconds, set pre and post-meal reminders to test, and more.

We are Exporters and Suppliers of Glucometer Accu chek Active at wholesale competitive prices from India.

Specifications :

FAST : Results in 5 seconds

MEMORY : 500 test results

Benefits and Features :-

- **Accurate:** Accurate¹ blood glucose measurement easy. No coding for even easier blood

glucose monitoring easy-to-read display (new typeface) intuitive handling with just 2 buttons.

- **Safe:** Underdose detection visual double check 10 seconds re-dose option
- **Clever:** Pre- and postmeal markers postprandial reminder 7,14,30 or now 90 day test averages USB interface for data transfer

- **Test principle :** Determination of glucose in fresh capillary blood by reflectance photometry.

Blood glucose values can be tested in whole blood or plasma.

The meter displays blood glucose values corresponding to those in plasma.

- **Measuring interval :** 10-600 mg/dL
- **Blood volume :** 1-2 unit
- **Measuring time :** Test strip in the meter: approx. 5 seconds
- **Test strip outside the meter:** approx. 08 seconds
- **Power supply :** 1 battery (type CR2032)
- **Battery life :** approx. 1000 tests or approx. 1 year
- **Automatic power-off :** after 30 or 90 seconds depending on operating status
- **Memory :** 500 test results with time and date, as well as 7, 14, 30 and 90 days averages
- **Temperature :**

1. During testing: +8 to +42°C
2. During storage without battery: -25 to +70°C
3. During storage with battery: -20 to +50°C

- **Atmospheric humidity :**

1. During testing: up to 85% relative humidity
2. During testing: up to 93% relative humidity

- **Altitude:** Sea level to 4000 m
- **Dimensions :** 97.8 x 46.8 x 19.1 mm
- **Weight :**

1. Without battery: approx. 46 g
2. With battery: approx. 50 g

- **Display :** 96-segment liquid crystal display (LCD)
 - **Interface :** USB (Micro B)
-

