

Labycare-Nebulizer Omron C 25S (SKSS-DP0067)



Nebulizer Omron C 25S Manufacturer by **LABYCARE** in India.

Description :

Labycare Nebulizer Omron C 25S is Compressor type Nebulizer

The compressed air from the compressor is forced through the nozzle and is mixed with the medication in the baffle and forced upwards. The medication is then turned into a fine aerosol and directed upwards by the baffle, where it is broken down into many smaller aerosol particles.

FOR EFFICIENT RESPIRATORY THERAPY

- Easy to operate
- High nebulization rate ensures optimized treatment duration
- Efficient delivery of medication
- Compact & travel convenient

Rating AC 230, 50Hz 1A

Power Consumption 75W

Nebulization Rate Approx 0.25ml/min (by weight loss)

Particle Size MMAD 4.21 um

MMAD = Mass Median Aerodynamic Diameter

Medication Tank Capacity 6ml maximum

Appropriate Medication Quantities 2ml minimum – 6ml maximum

Operating Temperature / Humidity +10°C to +40.0°C / 30% to 85% RH

Storage and Transport Temperature / Humidity / Air Pressure -20°C to +60°C / 10% to 95% RH / 700-1060hPA

Contents :

Compressor, Nebulizer Kit, Air Tube (180cm), Mouthpiece, Child Mask, 5 pcs Spare Air Filters, Carrying Case, Instruction Manual

We are Exporters and Suppliers of Nebulizer Omron C 25S at wholesale competitive prices from India.

Specifications :

Labycare Nebulizer Omron C 25S is Compressor type Nebulizer

The compressed air from the compressor is forced through the nozzle and is mixed with the medication in the baffle and forced upwards. The medication is then turned into a fine aerosol and directed upwards by the baffle, where it is broken down into many smaller aerosol particles.

FOR EFFICIENT RESPIRATORY THERAPY

- Easy to operate
- High nebulization rate ensures optimized treatment duration
- Efficient delivery of medication
- Compact & travel convenient

Rating AC 230, 50Hz 1A

Power Consumption 75W

Nebulization Rate Approx 0.25ml/min (by weight loss)

Particle Size MMAD 4.21 um

MMAD = Mass Median Aerodynamic Diameter

Medication Tank Capacity 6ml maximum

Appropriate Medication Quantities 2ml minimum – 6ml maximum

Operating Temperature / Humidity +10°C to +40.0°C / 30% to 85% RH

Storage and Transport Temperature / Humidity / Air Pressure -20°C to +60°C / 10% to 95% RH / 700-1060hPA

Contents :

Compressor, Nebulizer Kit, Air Tube (180cm), Mouthpiece, Child Mask, 5 pcs Spare Air

Filters, Carrying Case, Instruction Manual

We are Exporters and Suppliers of Nebulizer Omron C 25S at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com