

Labycare-Nebulizer Omron C 803 (SKSS-DP0069)



Nebulizer Omron C 803 Manufacturer by **LABYCARE** in India.

Description :

Labycare Nebulizer Omron C 803 is Compressor type Nebulizer

High nebulization rate

Small particle size- 3um

Compact size-ease of preparation, operation and storage

Low Noise

Rating (AC Adapter) : 100 -240V~50/60Hz 0.12 - 0.065A

Rating (Compressor Nebulizer) : 6V - 0.7A

Nebulization Rate : Approx 0.30ml/min (by weight loss)

Particle Size : MMAD approx 3um

Medication Tank Capacity 10 ml Maximum

Size : Medication Cup capacity 2 ml minimum - 10 ml maximum

Operating Temperature / Humidity : +10OC to + 40OC, 30% to 85% RH

Storage and Trasport Temperature / Humidity / Air Pressure : -20OC to + 60OC, 10% to 95% RH / 700 - 1060 hPa

Content :

Compressor, Nebulizer Kit, Air Tube (PVC, 100cm), Mouthpiece, Child Mask (PVC), 5 PC Spare Air Filters, AC Adapter, Instruction Manual.

We are Exporters and Suppliers of Omron Nebulizer at wholesale competitive prices from India.

Specifications :

Labycare Nebulizer Omron C 803 is Compressor type Nebulizer

High nebulization rate

Small particle size- 3um

Compact size-ease of preparation, operation and storage

Low Noise

Rating (AC Adapter) : 100 -240V~50/60Hz 0.12 - 0.065A

Rating (Compressor Nebulizer) : 6V - 0.7A

Nebulization Rate : Approx 0.30ml/min (by weight loss)

Particle Size : MMAD approx 3um

Medication Tank Capacity 10 ml Maximum

Size : Medication Cup capacity 2 ml minimum - 10 ml maximum

Operating Temperature / Humidity : +10OC to + 40OC, 30% to 85% RH

Storage and Trasport Temperature / Humidity / Air Pressure : -20OC to + 60OC, 10% to 95% RH / 700 - 1060 hPa

Content :

Compressor, Nebulizer Kit, Air Tube (PVC, 100cm), Mouthpiece, Child Mask (PVC), 5 PC Spare Air Filters, AC Adapter, Instruction Manual.

We are Exporters and Suppliers of Omron Nebulizer at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com