

Labycare-Otoscope Desco (SKSS-DP0073)



Otoscope Desco Manufacturer by **LABYCARE** in India.

Description :

Labycare Otoscope is also known as Auriscope. It is a medical device which is used to look into the ears and are used to screen for illness during regular check-ups and also to investigate ear symptoms.

An otoscope gives a view of the ear canal and tympanic membrane or eardrum. Otoscope Basic includes :

1 double lens otoscope head with 3 reusable speculas.

1 Battery handle with rheostat button.

We are Manufacturer, Exporters and Suppliers of Otoscope Desco at wholesale competitive prices from India.

Specifications :

- Labycare Otoscope is also known as Auriscope.
- Otoscope or Auriscope is a medical device which is used to look into the ears.
- Otoscopes are used to screen for illness during regular check-ups and also to investigate ear symptoms.
- An otoscope gives a view of the ear canal and tympanic membrane or eardrum Otoscope Basic includes :

1 double lens otoscope head with 3 reusable speculas.

1 Battery handle with rheostat button.

We are Manufacturer, Exporters and Suppliers of Otoscope Desco at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com