

Labycare-Vaginal Speculum Collin (SKSS00105)



Vaginal Speculum Collin _____

Description :

Labycare Vaginal Speculum Collin is made of stainless steel and is reusable.

Collin Vaginal Speculum is a gynecological instrument which is used for performing pelvic examinations, Pap Smears, and examining the cervix. It is a double blade duckbill speculum inserted with the blades closed. The blades are then opened by an external lever and screw.

There is no handle with the Collin device and the blades open laterally/horizontally versus vertically.

We are Manufacturer, Exporters and Suppliers of Vaginal Speculum Collin at wholesale competitive prices from India.

Specifications :

Labycare Vaginal Speculum Collin is made of stainless steel and is reusable.

Collin Vaginal Speculum is a gynecological instrument which is used for performing pelvic examinations, Pap Smears, and examining the cervix. It is a double blade duckbill speculum inserted with the blades closed. The blades are then opened by an external lever and screw.

There is no handle with the Collin device and the blades open laterally/horizontally versus vertically.

We are Manufacturer, Exporters and Suppliers of Vaginal Speculum Collin at wholesale competitive prices from India.