

Labycare-Vaginal Speculum Kristeller (SKSS00106)



Vaginal Speculum Kristeller ____

Description :

Labycare Vaginal Speculum Kristeller is made of stainless steel and is reuseable.

Kristeller Vaginal Speculum and Retractor Set are used in OB/GYN procedures.

They are commonly used for visual examinations because both are handheld, manual, non-locking instruments.

This set includes a 90mm x 36mm speculum and a 100mm x 32mm retractor. The sides of the speculum blade are curved upward slightly. The retractor has a smooth, flat blade. Both feature a slight hook at the bottom of the handle for pulling.

We are Manufacturer, Exporters and Suppliers of Vaginal Speculum Kristeller at wholesale competitive prices from India.

Specifications :

Labycare Vaginal Speculum Kristeller is made of stainless steel and is reuseable.

Kristeller Vaginal Speculum and Retractor Set are used in OB/GYN procedures.

They are commonly used for visual examinations because both are handheld, manual, non-locking instruments.

This set includes a 90mm x 36mm speculum and a 100mm x 32mm retractor. The sides of the speculum blade are curved upward slightly. The retractor has a smooth, flat blade. Both feature a slight hook at the bottom of the handle for pulling.

We are Manufacturer, Exporters and Suppliers of Vaginal Speculum Kristeller at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com