

## Labycare-Bacteriological Incubator Stainless Steel Chamber (SKSS362)



Bacteriological Incubator Stainless Steel Chamber Manufacturer by **Labycare** in India. It is used for storage of bacteria plate and bacterial culture growth at 37 degree Celsius. These incubators are fitted with heating temperature only. Therefore, these are also called heated incubators.

### **Description :**

Labycare Incubator has standard double wall fabrication. It is made out of stainless steel sheet and outer made out of thick mild steel sheet finished in white stoving enamel paint. The gap is filled with special glass wool for temperature insulation.

Double wall door has double window for observation in the chamber. Temperature ranges 5°C above ambient to 80°C + 2°C controlled by a thermostat.

L- shaped thermometer is provided for taking temperature readings. Panel is also provided with pilot lamps for mains & heaters with on-off switch and temperature Control knob.

Heating is achieved from 3 sides and perforated removable shelves are provided.

### **Specifications :**

- Double Walled Design
- Thick glass wool insulation to prevent heat loss
- Microprocessor based AUTO-TUNE PID Digital Temp Controller cum indicator.
- Temperature range : 50°C-250°C
- Temperature Accuracy : 0.5°C
- Temperature display : Digital LED Display
- Operates on : 230 volt, 50 Hz, Single Phase, AC Supply Special heavy duty stainless steel lock and door hinges.
- Heavy duty stainless steel trays for sample placement.
- Observation window to observe the sample without opening chamber door.
- Forced air circulation blower for Uniform Temperature Distribution Supplied complete with calibration certificate, IQ-OQ-PQ Document and test report.

- Calibration certificate with NPL Traceability for temperature sensor Digital controller.

- **Optional Accessories at Extra Cost :-**

I) Air Circulating fan

ii) Timer 0-24 hours

iii) Digital temperature Indicator Cum Controller

iv) Microprocessor PID digital temperature Indicator Cum Controller with timer

v) Digital Temperature Indicator

We are Manufacturer, Exporters and Suppliers of Bacteriological Incubator Stainless Steel Chamber at wholesale competitive prices from India.

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)