

## Labycare-Dry Bath Incubator (SKSS-336A)



Dry Bath Incubator Manufacturer by **Labycare in India**. It is a type of laboratory equipment that is used to heat samples. They are often used in molecular biology, microbiology, biochemistry and genetic applications. The capacity of these baths is measured in blocks.

### **Description :**

Labycare Dry Bath Incubator is also known as Dry Block Incubator.

Its Block has 25 Slots.

Incubator has micro controller design for precise temperature.

Incubator works on 220 Volt 50 Hz.

We are Manufacturer, Exporters and Suppliers of Dry Bath Incubator at wholesale competitive prices from India.

### **Specifications :**

- Model No : SKSS-336A
- Temperature Range : 5 to 75 degree Celsius
- Voltage : 220-240V, AC
- Frequency : 50 Hz
- Material : Mild Steel (Body)
- Automation Grade : Automatic
- Display Type : Digital

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)