

Labycare-Test Tube Vortex Shaker Manufacturer in India (SKSS-359A)



Test Tube Vortex Shaker Manufacturer by **Labycare** in India. It is a simple device used commonly in laboratories to mix small vials of liquid. In a vortex mixer, the motor drives the rubber cup in a circular motion to create a vortex, or spiral flow, in a sample. The idea is to place what needs to be mixed, such as a test tube, into the rubber piece so that it will rotate as well.

Description :

Test Tube Shaker is also known as Vortex Shaker. It is useful for accelerated mixing of solutions in test tubes. The housing is 150 x 150 mm. The speed is regulated by a speed control device. Short eccentric stroke motor imparts vigorous agitation to the tubes. The body of shaker is strong and is finished in stoving enamel paint. It is workable on 220V AC 50Hz single phase.

We are Manufacturer, Exporters and Suppliers of Test Tube Vortex Shaker at wholesale competitive prices from India.

Specifications :

Test Tube Shaker is also known as Vortex Shaker.

Vortex Shaker is useful for accelerated mixing of solutions in test tubes.

The housing is 150x150 mm.

The speed is regulated by a speed control device.

Short eccentric stroke motor imparts vigorous agitation to the tubes.

Body of shaker is strong and is finished in stoving enamel paint.

Vortex Shaker is workable on 220V AC 50Hz single phase.

We are Manufacturer, Exporters and Suppliers of Test Tube Vortex Shaker at wholesale competitive prices from India.

labycare

www.labycare.com | Email: labycare@gmail.com