

Labycare-Walk In Stability Chamber (SKSS365)



The Labycare Walk-In Stability Chamber is perfectly designed for the high requirements of stability studies and climatic tests. With stability rooms specially developed to meet FDA/ICH stability requirements generating exceptional control and uniformity of both temperature and humidity.

Description :

Walk In Stability Chamber powered by Labycare Technology for safe testing of specimens
European CE Certified Products Presafe DGM897.

Touch Screen and PLC based controller.

21 CFR Part 11 Ready Software.

Inbuilt USB and Data logger.

Stability testing for Pharmacy, Drugs, tablets, Cosmetic, Food, Beverages.

Specifications :

- Chamber Volume : 33cu ft 934 litres
- Chamber Dimension : 33w*27d*60h, 84*69*152 cm
- External Dimension : 41w*35d*90h, 104*89*229 cm
- Temperature : 2° C to 70° C
- Control Resolution : 0.1° C
- Uniformity at 37° C : $\pm 1^{\circ}$ C
- Heat-up Time : Amb -37°C : 7 Minutes
- Heat-up Time : 10°C -50°C : 22 Minutes
- Recovery time at 37°C : 1.25 Minutes
- Number of Shelves Supplied : 4
- Max Number of Shelves : 19

labycare

www.labycare.com | Email: labycare@gmail.com