

## Labycare-Non Destructive Moisture Meter (Imported) ()



Non Destructive Moisture Meter (Imported) Manufacturer by LABYCARE

### **Description :**

Non-Destructive Moisture Meter (Imported) is a high-precision device designed to measure moisture content in materials without damaging the sample. It uses advanced sensor technology to provide fast and reliable readings. Ideal for wood, grain, paper, and building materials, it ensures efficient quality control in laboratories, agriculture, and industrial applications.

### **Specifications :**

Commodity calibration	Vegetable seeds, Pulses, oil seeds, nuts, processed products, grass, for seeds etc.
Auto correction	Automatically corrects for entered test weights.
Display	LED display for moisture percentage. LCD display for commodity selection and user prompts.
Keypad	16 key Membrane keypad.
Auto Zero	The unit auto zeros itself when the cell is empty for short periods.
Sample weight	150-250gms.
Temp. correction	Automatic temperature correction.
Calibration	Calibrated to USDA basic oven standards.
Weight	15Kg approx.
Size	43cmx24cmx43cm (LxBxH) (approx.)
Power supply	220V, 50Hz.
Range	05% – 35%
Accuracy	±0.2%

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)