

Labycare-RICE POLISHER (SKSS-RLE2012)



RICE POLISHER Manufacturer by LABYCARE

Description :

The Rice Polisher is designed for the quick and efficient determination of rice polishing quality. It allows easy polishing of rice samples in a short time while maintaining high sample quality. The unit is self-contained and operates on the principle of abrasion and centrifugal force to polish the rice and separate the rice bran.

Construction: The Rice Polisher is designed for tabletop mounting and consists of an abrasive wheel and a mesh that polish the rice through abrasion while separating the bran. A collector cup is provided to collect the polished rice, while the rice bran is gathered in a separate collection drawer for further analysis. The entire unit operates using a single electric motor.

Specifications :

Capacity	100

Electric motor	0.5
Emery roller	30 m
Timer	Digi
Collecting pan	Mac
Step pulley	To s
Dimension	54 x

labycare

www.labycare.com | Email: labycare@gmail.com