

Labycare-Rice Analyzer Conveyor Based (SKSS-RLE2017)



Rice Analyzer Conveyor Based Manufacturer by LABYCARE

Description :

Perform rice sample quality control in just seconds with full traceability and objective results. Enhance your Quality Management System with technology that delivers accuracy, efficiency, and reliability.

- Automated Image Analysis: Inspects each grain individually to detect and classify damaged or defective kernels.

- Comprehensive Reporting: Generates detailed reports including data, charts, and annotated images of defective grains.
- Digital Traceability: Stores all sample results for future reference, ensuring complete recordkeeping.
- Boosts Lab Productivity: Increases accuracy, reduces human error, and speeds up workflows in quality control labs.

Specifications :

Multispectral analysis RGB	Infrared 850 / 940nm.
Analysis and classification	Defects and foreign bodies.
2D Granulometry	Length and width measurement.
Accurate reports	with data, graphics and images.
Speed	60 sec. per sample of 100 grs.
Results tickets	Instant ticket emission including main data.

labycare

www.labycare.com | Email: labycare@gmail.com