

Labycare-PALM DEHUSKER WITH BROKEN SEPARATOR (SKSS-RLE2019)



PALM DEHUSKER WITH BROKEN SEPARATOR Manufacturer by LABYCARE

Description :

The Palm Dehusker with Broken Separator is a compact and reliable device designed for manually dehusking small paddy samples. Constructed from high-quality ABS and rubber, it offers durability and precision in operation. The dehusker features grooved rubber surfaces on both the upper and lower sections, ensuring efficient husk removal. With an integrated broken separator, the unit comes with three plastic screens featuring 1.6 mm, 3.2 mm, and 4.8 mm indents for effective separation of brokens from head rice. This makes it a valuable tool for laboratories, research facilities, and quality control processes, providing accurate sample analysis and grading.

Specifications :

- Designed for manual dehusking of small paddy sample.
- Made of special quality ABS and rubber.
- Upper part and one lower part, both fitted with grooved rubber.
- The broken separator is supplied with three plastic screens having 1.6 mm, 3.2 mm, and 4.8 mm indent to separate brokens from head rice.

labycare

www.labycare.com | Email: labycare@gmail.com