

Labycare-ANALOG GRAIN VERNIER (SKSS-RLE2021)



ANALOG GRAIN VERNIER Manufacturer by LABYCARE

Description :

Analog Grain Vernier (Dial Type) is a precision mechanical instrument designed to measure rice grain thickness accurately. It features a clear circular dial indicator and stable grain holding slot for consistent readings. Ideal for rice mills, laboratories, and quality control units requiring reliable manual measurement.

Specifications :

Range	0 to
Accuracy	0.01
Dimensions	145
Weight	345

labycare

www.labycare.com | Email: labycare@gmail.com