

## Labycare-ILLUMINATED PURITY WORK BOARD (SKSS-RLE2022)



ILLUMINATED PURITY WORK BOARD Manufacturer by LABYCARE

### **Description :**

The Illuminated Purity Works Board is designed for the examination of grains, broken rice, and other similar samples. A light placed beneath the glass surface shines through, enabling a detailed and critical inspection of the sample. This equipment is primarily used by grain inspectors and seed analysts. The illuminated light source creates a well-lit work area on the smooth glass surface, ideal for grain testing.

It is equipped with a magnifying glass that covers the entire surface of the opal glass. The large magnifying lens has a diameter of 100 mm (4 inches) and a focal length of 20 cm. The board features an on/off switch and two 25W candle-type bulbs. The opal glass also includes a scale for measuring the length of the grains.

<b>Specifications :</b>	Testing area opal glass size	11?
	Magnifying glass	20 d
	Bulb	2 no
	Drawer	1 no

---

labycare

[www.labycare.com](http://www.labycare.com) | Email: [labycare@gmail.com](mailto:labycare@gmail.com)