

Labycare-FLEXIBLE ARM ILLUMINATED MAGNIFIER (SKSS-RLE2023)



Flexible Arm Illuminated Magnifier Manufacturer by LABYCARE

Description :

See every detail clearly with our Flexible Arm Illuminated Magnifier. Designed for hobbyists, professionals, and anyone needing sharper vision, it combines high-quality magnification and adjustable LED lighting to make detailed work easy and enjoyable.

<u>Specifications :</u> Illumination ?????	A 22 watt circular fluorescent tube housed in a tough ABS shed around the lens provides light that is shadow less, glare free, cool and of daylight colour.
Versatile stand ?????	Easy positioning of lens / light at any point in any orientation is made of possible by a well balanced stand with feather touch movement that has a reach of 900 mm from the base.

Documents : <images/babfaed7e5eafaafa3e21ae78d56c4d.jpg>

labycare

www.labycare.com | Email: labycare@gmail.com